

“ Simple straight-forward flow testing system ”

The FT-Flowtester is a simple and straight forward system to test the difference between the flow before and after coating and/or drying/conditioning.



After placing seed in the hopper you can set the opening time of the valve. After pressing the green button the valve opens during the set time. You can put a holding bucket or container underneath the slide. This platform can be easily change in height, so several different containers can be used.

When you separately weigh the dosed amount in the container/bucket you can see if the flow characteristics are changed.

Product characteristics:

- Hopper for approx. 30l of seeds
- Pneumatic valve for dosing
- Timer to set opening duration
- Connections : 1~110-230+N Vac 50hz
: Compressed air for the valve (6 bar)

FT-Flowtester

